

	Table of Contents
--	-------------------

Part S1	Hardware description	S2
	Description of main board	S3
	Description of key board	S5
	Description of spo2 board	S6
Part S2	Performance verification	S7
	Required Equipment	S7
	Pulse Oximetry Functional Test items	S7
	Pulse Oximetry Functional Tests	S7
Part S3	Troubleshooting	S9
	Password for maintenance menu	S9
	Troubleshooting	S9
Part S4	process for software up-grade	S11
	PC to device	S11
	To prepare download program(PC)	S12
	Act & Set up FlashMagic	S12
	Executing upgrade of program	S15
	Double checking update status by	S17
Part S5	Board circuit	S18

Part S1 Hardware Description




Fig. Whole Block Diagram

This product consists of three boards including Main, Key and Analog (SpO2) boards. On individual board, Main board has general system control and power supply parts and Key board is equipped with function buttons and a sound amplifier to generate sound. Lastly Analog (SpO2) board can detect analog signals.

Description of Main Board




Fig. Block Diagram of Main board

■ CPU

● Structure of CPU

- 16/32-bit ARM7TDMI-S Micro Controller
- 16 kB One Chip RAM
- 128/256 kB One Chip Flash Program Memory
- 10-bit A/D Converter
- Two 32-bit Timer, Real-time Clock and Watchdog

● Functions of CPU

- To control brightness of LCD.
- To display saturation level of oxygen and pulse frequency for the patient.
- To store values measured from the patient.
- To check button pressed and execute the functions upon it.

- To calculate saturation level of oxygen and pulse frequency.
- To generate alarm, pulsation sound and button pressing sound.
- To check remaining battery capacity and control power supply of the product.

■ Charging Circuit

- To check voltage and current of Li-ion battery (4.2V) in real-time and charge it by applying Voltage.
- To check voltage, current and charging time constantly to ensure safe charging.

■ Power Supply Circuit

- If input voltage goes down under 5V, a step-up converter operates to generate stable voltage +5V so that alkaline battery can be used.
- To generate +3.3V and +1.8V to be used on Main board.

■ SRAM

- It consists of 1MByte RAM, which stores values measured from the patient and settings for the product.

■ Real-time Clock (RTC), Watchdog, SRAM Controller

- bq4802 includes all functions for Real-time Clock (RTC), Watchdog, SRAM Controller, etc.
- Real-time Clock (RTC) is used to set the current time.
- Watchdog carries out reset when CPU doesn't work so that it can operate again.
- SRAM controller protects SRAM from storing any abnormal data when turning on/off power supply at the beginning or the power is unstable.

Description on Key Board





Fig. Block Diagram of Key Board

- Power Control Logic sends signals to power circuit and CPU on Main Board to turn on/off the product when pressing power button.
- It receives power applied from Main board and generates +3.3V to allow power control logic to operate.
- It has six buttons except power button, of which pressing status are transferred to CPU.

Description on Analog (SpO2) Board



Block Diagram of Analog (SpO2) Board

■ CPU (AVR)


- It consists of 8bit micro controller to carry out functions such as communication with Main CPU, current control, signal separation, A/D converter control, etc.

■ Procedures to Extract Signals

- To control current driver with value set by CPU to apply it to the sensor.
- When a finger enters into the probe (sensor), it amplifies the signal and extracts the signal through low pass and high pass filters. The extracted signals are divided into RED and IR signals through MUX.
- The divided signals are fitted with the protocols set by CPU, and sent to CPU on Main Board through serial port.

Part S2 Performance verification

-. Required Equipment

Equipment	Description
Pulse Oximetry Sensor	Charmcare reusable Sensor
Pulse Oximetry Functional Tester	Fluke Index2 


**Note**

A functional tester cannot be used to assess the accuracy of a pulse oximeter monitor. However, it can be used to demonstrate that a particular pulse oximeter monitor reproduces a calibration curve that has been independently demonstrated to fulfill a particular accuracy specification.

-. Pulse Oximetry Functional Test items

Item	Description
Pulse rate(BPM)	The test procedure simulates an sensor attached to a patient indicating 60 BPM and 200 BPM.
SpO2	The test procedure simulates an sensor attached to a patient indicating 66% blood oxygen saturation and 96% blood oxygen saturation.

-. Pulse Oximetry Functional Tests

1	With the monitor turned off, connect the pulse oximetry sensor to the sensor port.
2	The Index2 tester set 96% blood oxygen saturation and 60 BPM (Main Menu1 → SIM → MAN → SPO2 and RATE)
3	Connect the Index2 tester to the other end of the sensor
4	Turn on the monitor by pressing the  Power button.

5	Check SpO2 and pulse rate measurement values. -.%SpO2 indication between 94 and 98 inclusive. -.BPM indication between 57 and 63 inclusive.
6	The Index2 tester set 66% blood oxygen saturation and 200 BPM (Main Menu1 → SIM → MAN → SPO2 and RATE)
7	Check SpO2 and pulse rate measurement values -.%SpO2 indication between 63 and 69 inclusive. -.BPM indication between 197 and 203 inclusive.

Part S3 Troubleshooting

-. Password for maintenance menu

Password is required to access maintenance. The password is **1289**.




-. Troubleshooting


Symptom	Problem	Action
Power does not turn on.	In case a cable between main board and key board is unstable or broken.	Check connection status of cable and its own status. (J1 on main board and J1 on key board)
	In case battery connection has any problem.	In case of using disposable battery, check the cable connected to the main board.
		In case of using rechargeable battery, check charging status or the cable connected to the main board.
Power turns on but other button does not work.	In case a cable between main board and key board is unstable or broken.	Check connection status of cable and its own status. (J3 on main board and J1 on key board)
Screen does not appear.	No initial setting	▲Press and hold (upper) button and ▼(lower) button, and turn on power.
	In case a cable between LCD and main board is unstable.	Check connection status of cable. (J3 on main board)
No sound.	In case connection between speaker and key board is broken or unstable.	Check connection status of cable and its own status. (J2 on key board)
	In case a cable between main board and key board is unstable or broken.	Check connection status of cable and its own status. (J1 on main board and J1 on key board)
As screen did not appear, initial setup was carried out to display the screen. However, when turning power again,	Connect on backup battery has any problem or the battery has been exhaustively discharged.	Check or replace backup battery. (BT1 on main board)

the screen does not appear again.		
Screen displays [SPO2 FAULT].	Main board and analog (SPO2) board cannot communicate each other.	Check connection status of main board and analog (SPO2) board. (J8 on main board J5 on analog board)
Screen displays [CLOCK ERROR].	Current time setting exceeds allowable range.	Set up the current time again on main board.


Part S4 Process for software up-grade

PC to device


Picture	Description
	<p>prepare the download cable – RS232 (Serial 9pin).</p>
	<p>Check the existence of serial port(COMx) in the back of your computer.</p>
	<p>Connect the download cable to serial port(COMx).</p>

**To prepare download program (PC)**

1. Install the file of "FlashMagic_VerX.X.exe" Charmcare offered..

**Act & Set up FlashMagic**


1. Do the following instructions to set up as the following picture.



[Pic. 1]

Step1) Communication

- 1) Press the button of "Select Device" → Select ARM7 → Select "LPC2214"



[Pic. 2]

2) Set up COM port.

COM Port : COM1
 Baud Rate :115200
 Interface :None(ISP)
 Oscillator(MHz) :12

Step2) Erase

Set up as [Pic. 1].

Step3) Hex File

Press the button of "Browse", and select "**0A4MAINARm1xxx.hex**" file Charmcare offered.


Step4) Options

Set up as [Pic. 1].


Executing update of program

- 1) Check whether the cable connects normally PC and the device.
- 2) Power on the device.
Notice) The screen of the device is not powered on at this moment.
- 3) Execute the program, "FlashMagic.exe", and check the set up status.
- 4) Press Step5) "**Start**" button of "FlashMagic.exe" for updating.

Its progress status appears at the bottom of screen as the following picture.



[Pic. 3]




[Pic. 4]

- 5) If Update is done, the program status disappears.
- 6) Unplug the cable connected to the device, and then press the power button roughly for 2 seconds to power off.





Double checking update status by

- 1) Check whether the cable is connected to the device. If it's connected, unplug it.
- 2) Power on by pressing the power button of the device.
- 3) Press the Menu button to move to the main Menu screen.
- 4) Version No. and the date for this programming appears at the bottom of Menu screen as the following picture.




Part S5 Board circuit

No	Board List
1	Main board
2	Spo2 board
4	Key board



- TP10 TEST POINT ○ SW1
- TP11 TEST POINT ○ SW2
- TP12 TEST POINT ○ SW3
- TP13 TEST POINT ○ SW4
- TP14 TEST POINT ○ SW5
- TP15 TEST POINT ○ SW6
- TP16 TEST POINT ○ POWER_EN/ALBE
- TP17 TEST POINT ○ B_POWER_SW
- TP18 TEST POINT ○ POWER_ON/OFF
- TP19 TEST POINT ○ POWER_EN/ALBE
- TP20 TEST POINT ○ ADAPTOR_INPUT
- TP21 TEST POINT ○ +4.5.0DV
- TP22 TEST POINT ○ PWM_SOUND
- TP23 TEST POINT ○ +3.0DV



Spo2 board

