


# CLINICAL SPOTLIGHT


## Table-Mounted Retraction Systems for Urology

Mediflex<sup>®</sup>, a trusted brand for over 50 years, offers retraction systems to fit specific patient size or procedural needs. A variety of table mounted options are available which can be used for other open, minimally invasive or endoscopic procedures.



### Mini-Bookler<sup>®</sup> TMR System

- Miniature components are less traumatic and offer a lower profile
- Split ring allows for frame adjustment according to the surgical site size and location
- Frame aids in supporting/retracting adjacent anatomy for maximum field exposure
- Frame notches permit use of disposable Surgical Stays for superficial retraction
- StrongArm<sup>™</sup> or FlexArm<sup>™</sup> enables versatile, low-profile retractor positioning



Bulbar Urethroplastic Surgery

**Assoc. Prof. Mazhar Ortac, MD, PhD**

Istanbul University Capa Med School - Urology Dept.



Bulbar Urethroplastic Surgery

**Assoc. Prof. Dr. Fatih Yanaral – Urologic Surgeon**

Haseki State Hospital - Istanbul, Turkey

### Bookler<sup>®</sup> TMR System

- Mediflex<sup>®</sup> Bookler<sup>®</sup> components are manufactured of the same materials and are dimensionally equal to leading systems – making component replacement seamless with guaranteed interchangeability
- Wide array of blade configurations
- Protected edge Malleable Blades
- CEU accredited inservice for inspection and processing available

**For more information, contact us at:**  
631.582.6424 or [info@mediflex.com](mailto:info@mediflex.com)

**[www.mediflex.com](http://www.mediflex.com)**